

References:

1. Renstrom A, Larsson PH, Malmberg P, Bayard C. A new amplified monoclonal rat allergen assay used for evaluation of ventilation improvements in animal rooms. J Allergy Clin Immunol. 1997 Nov;100(5):649-55.
2. Renstrom A, Gordon S, Larsson PH, Tee RD, Newman Taylor AJ, Malmberg P. Comparison of a radioallergosorbent (RAST) inhibition method and a monoclonal enzyme linked immunosorbent assay (ELISA) for aeroallergen measurement. Clin Exp Allergy. 1997 Nov;27(11):1314-21.




Rat n 1 ELISA kit (RUP-6/RUP-1)

Product Code: EL-RN1

Lot Number: XXXXX

Sample Curve:



www.inbio.com

Indoor Biotechnologies, Inc.
700 Harris Street,
Charlottesville, VA 22903
United States

Indoor Biotechnologies Ltd.
Vision Court
Caxton Place
Cardiff, Wales CF23 8HA
United Kingdom

Indoor Biotechnologies,
India Private Limited
Bangalore Bioinnovation Centre,
BioTech Park, Electronic City Phase 1
Bangalore-560100, India

Tel: (434) 984-2304
Fax: (434) 984-2709
Email: mail@inbio.com

Tel: +44 (0) 29 2167 4640
Email: Info@indoorbiotech.co.uk

Tel: +91-9901722009
Email: info@inbioindia.com

Content:

Vial 1 (red top) 100 μ L
Monoclonal antibody RUP-6

Vial 2 (white top) 400 μ L
Universal Allergen Standard
Concentration: 1000ng/ml Rat n 1

Vial 3 (brown) 100 μ L
Biotinylated monoclonal antibody RUP-1
Dilute: 1:1000 for use

Storage: The ELISA kit should be stored at 4°C

For research and commercial use in vitro: not for human in vivo or therapeutic use.

An InBio™ product

Certificate of Analysis

Monoclonal Antibody: RUP-6
Immunogen: Rat n 1
Isotype: Mouse IgG1
Specificity: Binds to an epitope on rat *Rattus norvegicus* urinary allergen, Rat n 1.
Purification: Produced in cell culture and purified by affinity chromatography using Protein G.
Concentration: 1.0mg/ml in phosphate buffered saline. Based on A280 for IgG (1.42= 1mg/ml) 0.22µm filtered, 0.02% sodium azide.
Lot Number: xxxxx

Monoclonal Antibody: RUP-1
Immunogen: Rat n 1
Isotype: Mouse IgG1
Specificity: Binds to an epitope on rat *Rattus norvegicus* urinary allergen, Rat n 1.
Purification: Produced in cell culture and purified by affinity chromatography using Protein G.
Concentration: 1.0mg/ml in phosphate buffered saline. Based on A280 for IgG (1.42= 1mg/ml) 0.22µm filtered, 0.02% sodium azide.
Biotinylation: Biotinylated and titrated for use in ELISA at 1/1000 dilution.
Lot Number: xxxxx

Allergen Standard: Universal Allergen Standard
Composition: A formulation of eight purified natural allergens prepared in 1% BSA/50% glycerol/PBS, pH 7.4
Lot Number: xxxxx

Concentration:

Universal Allergen Standard	Protein Measurement	Concentration (ng/ml)
Der p 1	Amino-acid analysis	2500
Der f 1	Amino-acid analysis	2500
Der p 2	Amino-acid analysis	1000
Fel d 1	Amino-acid analysis	1000
Can f 1	Amino-acid analysis	2500
Rat n 1	Amino-acid analysis	1000
Mus m 1	Amino-acid analysis	250
Bla g 2	Amino-acid analysis	2500

ELISA Protocol for Rat n 1.

- Coat polystyrene microtiter plates (NUNC Maxisorp Cert. NUNC catalog # 439454) with 100µl mAb RUP-6 at 10µl/10ml, i.e. 1/1000 dilution of stock, in 50mM carbonate-bicarbonate buffer, pH 9.6, incubate overnight at 4°C.
- Wash wells 3x with PBS-0.05% Tween 20, pH 7.4 (PBS-T). Incubate for 30 min. at room temperature with 100µl/well of 1% BSA, PBS-T. Wash 3x with PBS-T.
- Use doubling dilutions of the Universal Allergen Standard to make a control curve ranging from 100 - 0.2ng/ml Rat n 1: Pipette 20µl UAS into 180µl 1% BSA, PBS-T into wells A1 and B1 on the ELISA plate. Mix well and transfer 100µl across the plate into 100µl 1% BSA, PBS-T diluent to make 10 serial doubling dilutions. Wells A11, B11 and A12, B12 should contain only 1% BSA, PBS-T as blanks.
- Add 100µl of diluted allergen samples and incubate for 1 hour at room temperature. House dust extracts for Rat n 1 analysis are routinely diluted two-fold from 1/10-1/80. Other sample types, like air filter extracts and allergen extracts, may require different dilutions.
- Wash wells 3x with PBS-T and add 100µl diluted biotinylated anti-Rat n 1 mAb RUP-1. The antibody solution should be diluted 1/1000 in 1% BSA, PBS-T. Incubate for 1 hour at room temperature.
- Wash wells 3x with PBS-T and add 100µl diluted Streptavidin-Peroxidase (Sigma S5512, 0.25mg reconstituted in 1ml distilled water). The reconstituted Streptavidin should be diluted 1/1000 in 1% BSA, PBS-T. Incubate for 30 minutes at room temperature.
- Wash wells 3x and develop the assays by adding 100µl 1mM ABTS in 70mM citrate phosphate buffer, pH 4.2 and 1/1000 dilution of H₂O₂. Read the plate when the absorbance at 405nm reaches 2.0-2.4.

Notes:

Universal Allergen Standard is recommended for immunoassay calibration purposes only. Not recommended for in-vitro antibody measurements, T cell studies, immunization purposes, or other uses.

Buffer recipes, storage conditions and a list of frequently asked questions can be found under "Protocols" on our web site: www.inbio.com.

For research and commercial use in vitro: not for human in vivo or therapeutic use.